
**Direct Reading
Echelle ICP**

**High Throughput Sample
Analysis with the Leeman
Labs Direct Reading Echelle
ICP**

Application Note 1019

High Throughput Sample Analysis with the Leeman Labs Direct Reading Echelle ICP

Authored by:
Manuel Almeida, Ph.D.

Leeman Labs, Inc.
6 Wentworth Drive
Hudson, NH 03051
U.S.A.
Tel: 603-886-8400 • Fax: 603-886-9141
Website: <http://www.LeemanLabs.com>

Leeman Labs (UK) Limited
University of Warwick Science Park
Barclays Venture Centre
Sir Williams Lyons Road
Coventry, CV4 7EZ
United Kingdom
Tel: +44-1203-690490
Fax: +44-1203-690404

Leeman Labs Germany GmbH
Zum Hammerseifen 46
D-57223 Kreuztal
Germany
Tel: +49-2732-4020
Fax: +49-2732-4010

Leeman Labs Benelux B.V.
van Aerschotstraat 56
4761 GE Zevenbergen
The Netherlands
Tel: +31-168-370241
Fax: +31-168-370242



Value Through Innovation

Introduction

Many high throughput labs are required to analyze hundreds of samples a day. Included in this category are oil, geological, agricultural, food and environmental labs. If the element suite is large enough, then the only ICP capable of handling the sample load is a fixed-channel simultaneous. A Simultaneous Leeman Labs Direct Reading Echelle (DRE) ICP combines the high resolution and dispersion of echelle optics with high gain detectors to provide the fastest sample throughput possible. ICPs using array type detectors simply cannot match the speed of the PMT based DRE, particularly when the samples are at low concentrations in a complicated matrix. Low UV sensitivity combined with the lack of signal multiplication works against the array-based systems in these high throughput applications.

The sample analysis time of a simultaneous (multi-channel) DRE is approximately one minute, including uptake and rinse. The Direct Reading Echelle (DRE) High-Throughput (HT) ICP with a new version of software, some changes to the basic analysis

sequence and modifications to sample introduction hardware, has reduced analysis time by almost 50%. This results in a dramatic improvement in productivity for the high volume lab — running nearly twice as many samples as they currently do in a given day.

Instrumentation

In order to determine the increase in sample throughput, the modifications were made to a multi-channel DRE operating in a high volume oil lab. Performing wear metals in oil analysis, the lab analyzes about 300 samples per day. Each sample is analyzed for 11 elements, which are listed in Table 1 along with wavelengths and background correction points. The integration time used was 3 seconds with a 1 second background correction. Without modifications, the analysis time for a single sample was 72 seconds.

The high throughput modifications were tested with real samples submitted to the lab for analysis. The test was conducted by first running the day's

| Element | Wavelength, nm | Background Correction Wavelength, nm |
|---------|----------------|--------------------------------------|
| Cu | 324.754 | 324.717 |
| Fe | 259.940 | 259.910 |
| Cr | 267.716 | 267.685 |
| Al | 308.215 | 308.250 |
| Si | 288.158 | 288.125 |
| P | 220.353 | 220.378 |
| Mo | 277.540 | 277.508 |
| Na | 589.592 | 589.660 |
| Sn | 189.926 | 189.904 |
| Ni | 231.604 | 231.577 |
| Zn | 206.200 | 206.224 |

Table 1 Elements, Wavelengths and Background Correction Points

| Table 2 Analysis Conditions | |
|--------------------------------|------------|
| RF Power | 1.1 kW |
| Coolant Flow | 18 LPM |
| Auxiliary Flow | 0.5 LPM |
| Nebulizer Pressure | 45 PSI |
| Pump Rate | 1.1 ml/min |

sample load using the modified sample introduction system but none of the software changes. The sample load consisted of 300 samples plus periodic QC standards, which made the total number of analyses to be 320. The analytical conditions used are listed in Table 2. The entire sample load was run without a single QC failure, indicating the modified sample introduction system and the V-groove nebulizer had the necessary stability.

Once the lab's work was complete, the software modifications were made. These consisted of routines designed to accelerate the uptake and rinse cycle of each sample. By moving the sample tip back to the rinse before analysis is complete and increasing the speed of the peristaltic pump, 30 seconds could be removed from the analysis time. With all the modifications, the analysis time was reduced to 42 seconds per sample.

Results

The first batch of samples analyzed by the lab was re-analyzed with the high throughput modifications. Comparisons of some of the results obtained are shown in Table 3. These data are corrected for the 1:4 dilution performed on the original oil samples. In addition, these data are rounded to the nearest integer, being the standard practice for this laboratory.



The data indicates there is no difference in accuracy between the standard and

high throughput set-ups. The data for all elements correlated extremely well between the two analysis runs. Figures 1 and 2 demonstrate the excellent correlation between the standard and high throughput mode for Fe ($r = 0.9997$) and Pb ($r = 0.9986$), respectively. All the elements demonstrated similar degrees of correlation over the entire concentration range found in the samples. The accelerated uptake and rinse times are sufficient to prevent any carryover between samples. This is extremely important since real samples such as these may differ significantly in concentration. During the second set of analyses, the same QC standard frequency was used and no failures occurred. This indicates the modified set-up is just as stable as the standard one. And best of all, the analysis time was reduced to 42 seconds a sample. At this analysis rate, a total of approximately 550 analyses could be completed in the time taken to run the 320 under the original set-up.

| Standard and High Throughput Mode Results Comparison | | | | | | | | | |
|--|------|------|------|-----|------|------|------|------|---------|
| Sample | 1002 | | 1004 | | 1006 | | 1008 | | Element |
| | STD | HT | STD | HT | STD | HT | STD | HT | |
| Cu1 | 15 | 13 | 16 | 16 | 3 | 3 | 2 | 2 | |
| Fe2 | 4 | 4 | 33 | 33 | 10 | 9 | 17 | 17 | |
| Cr4 | 0 | 0 | 2 | 2 | 0 | 0 | 0 | 0 | |
| Al3 | 2 | 1 | 14 | 12 | 3 | 3 | 0 | 0 | |
| Si2 | 6 | 5 | 26 | 25 | 7 | 7 | 3 | 3 | |
| Pb1 | 1 | 1 | 1 | 1 | 2 | 2 | 1 | 1 | |
| Mo3 | 1 | 2 | 2 | 2 | 3 | 3 | 4 | 4 | |
| Na2 | 4 | 2 | 1 | 0 | 3 | 2 | 2 | 2 | |
| Sn1 | 1 | 1 | 1 | 1 | 1 | 1 | 1 | 1 | |
| Ni3 | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 0 | |
| Zn3 | 1108 | 1084 | 915 | 896 | 1046 | 1048 | 1295 | 1295 | |
| Sample | 1027 | | 1029 | | 1031 | | 1033 | | Element |
| Element | STD | HT | STD | HT | STD | HT | STD | HT | |
| Cu1 | 4 | 4 | 3 | 3 | 0 | 0 | 7 | 7 | |
| Fe2 | 67 | 66 | 67 | 65 | 3 | 3 | 15 | 16 | |
| Cr4 | 1 | 1 | 1 | 1 | 0 | 0 | 1 | 1 | |
| Al3 | 2 | 3 | 2 | 2 | 1 | 1 | 7 | 8 | |
| Si2 | 16 | 18 | 16 | 17 | 11 | 12 | 18 | 16 | |
| Pb1 | 1 | 1 | 1 | 1 | 0 | 0 | 1 | 1 | |
| Mo3 | 1 | 2 | 1 | 2 | 1 | 1 | 98 | 97 | |
| Na2 | 0 | 0 | 0 | 0 | 0 | 0 | 2 | 2 | |
| Sn1 | 0 | 1 | 1 | 1 | 0 | 1 | 1 | 1 | |

Table 3 Standard and High Throughput Data Comparison

STD = Standard Analysis Speed
HT = High Throughput Mode



Value Through Innovation

**LEEMAN
LABS, INC**

6 Wentworth Drive
Hudson, NH 03051
Tel.: 603-886-8400
Fax.: 603-886-9141

