

PRESS RELEASE

March 2006

Teledyne Leeman Labs introduces a new transportable simultaneous ICP, the Prodigy - Mobile Lab

Hudson, NH - Teledyne Leeman Labs, a leading manufacturer of analytical instrumentation for elemental analysis, announces the release of a new transportable simultaneous ICP spectrometer, the Prodigy-Mobile Lab. By providing accessibility to laboratory quality ICP-AES instrumentation where you need it, the Prodigy-Mobile Lab can accelerate the collection of data, reduce the cost of producing it and help you make more informed decisions about site characterization.

Our engineers worked closely with truck and trailer manufacturers as well as remediation companies to determine the specific hazards associated with on-site transportable operation. Prodigy-Mobile Labs' integral suspension system works in conjunction with parent vehicle's suspension, to isolate the instrument from drive train, road and shock hazards. The result is laboratory quality data for in-field operation.

For information on the Prodigy-Mobile Lab, contact Teledyne Leeman Labs, 6 Wentworth Drive, Hudson, NH 03051. Telephone (603) 886-8400, Fax: (603) 886-9141 or email: LeemanLabsinfo@Teledyne.com.

For more information on our other products or services visit our website at: www.LeemanLabs.com.